
## COVID CARE


Reimagining the health checkup procedures in post covid era

# Our inspiration for this project

As the world will reopen after the pandemic, there will be a necessity to follow the covid SOP in public places and the need for a one-stop solution for checking all the vitals will surge. Our project COVID CARE was made to serve this need. It focuses on providing maximum services to track the health of a person and maintain a database of the health of people checking in a public place.



## COVID CARE: What it does




Provides an instant report of a person's pulse rate, oxygen levels and body temperature.



Access your health records at a public check-in on mobile app.




Maintains a database of the health of all the people checking in at a public place.



Sends alert in case of adverse changes in health parameters




The user of the mobile app can also add health records.



#### Sensors used to detect vitals in Proteus




Pulse sensor



Temperature Sensor (Tmp36)

Note: Here we have used a pulse sensor instead of a pulse oximeter sensor, because library of the pulse oximeter sensor for proteus and other softwares was not available.

# Future Aspects of Hardware Implementation



In real life scenario we will be using MAX 30100 pulse oximeter sensor module. It's an optical sensor that derives its readings from emitting two wavelengths of light from two LEDs - a red and an infrared one - then measuring the absorbance of pulsing blood through a photodetector. This particular LED color combination is optimized for reading the data through the tip of one's finger. The signal is processed by a low-noise analog signal processing unit and communicated to the target MCU through the mikroBUS I2C interface.

# What's next for COVID CARE

- Ol UI/UX designs of mobile app with accessibility features for specially abled
- Integration with fitness apps and devices like Fitbit
- Adding Wi-fi and Bluetooth functionality in device to connect with mobile app
- Using pulse oximeter sensor instead of pulse sensor in hardware part

### Our team

Shivakshi Sheel Srivastava



Priyanshi Chandra



Nishi Sharma



Nisha Prakash



## Thank you